

XXL Prosthetic Care – Heavy-Duty Line: Components for users with body weight of up to max. 225 kg



XXL Prosthetic Care: Components for users with body weight of up to max. 225 kg

Streifeneder ortho production GmbH's heavy duty product line covers components for maximum body weight of up to 225 kg (496 lbs) and has been developed to restore the mobility of heavy weight users. With these components, the orthopaedic technician is able to build challenging prosthetic devices and at the same time, can depend on consistently excellent quality.

Features

- prosthetic care for users with a body weight of up to 225 kg (496 lbs)
- low construction weight plus maximum rigidity due to high quality titanium and aluminium
- consistently excellent quality "made by Streifeneder ortho.production"
- approved according to ISO 10328 standards



Screw Adapter -Titanium-, Ø 34 mm 11A250/T

- eff. height: 16 mm
- for tubes with Ø 34 mm
- adjustment core adapter suitable for standard modular system
- with 4 adjusting screws 22A2/M8x12mm, zinc-plated
- material: titanium; weight: 131 g; total height 55 mm
- Modular system up to max. 225 kg – with straight, non-adjustable construction!



Screw Adapter -Titanium-, Ø 34 mm 11A225/T

- eff. height: 21 mm
- tube clamp Ø 34 mm
- adjustment core adapter exclusively suitable for S.A.C.H. foot adapter 16A7/M14 and socket adapter 10A7/T
- with 4 adjustment set screws 22A2/M10x16
- material: titanium; weight: 167 g; overall height 55 mm
- adjustment possible in connection with 225 kg pylons!



© starush / fotolia



Socket Adapter 10A7/T

- with adjustment core; overall height 52 mm
- for BK socket
- extra wide notch
- for direct lamination
- adjustment core exclusively suitable for screw adapter 11A225/T
- material: titanium; weight: 213 g; eff. height: 23 mm



Light Alloy Tube 28A3/A300

- light alloy tube \varnothing 34 x 6 mm
- silver anodized
- suitable for adapters with tube holder \varnothing 34 mm
- material: aluminium; weight: 460 g; length: 300 mm



S.A.C.H. Foot Adapter -Titanium- 16A7/M14

- complete with screwing
- steel screw with fine thread M14x1,5
- use only in combination with screw adapter 11A225/T



S.A.C.H. Foot 30A225/..

- with fixed S.A.C.H. Foot Adapter 16A7/M14
- height incl. S.A.C.H. Foot Adapter: 80 – 90 mm
- heel height: 10 mm \pm 2 mm
- available in sizes 25 – 30
- activity level (walky): 1

The component with the lowest maximum (body) weight in the construction of a prostheses is decisive for the maximum load of the entire prosthesis.

Streifeneder
ortho.production GmbH
Moosfeldstraße 10
82275 Emmering
Germany
T +49 8141 6106-0
F +49 8141 6106-70
export@streifeneder.de
www.streifeneder.com/op



300W104/E 2017-03