

OmoSupport.easy Shoulder Joint Sling Orthosis for immobilisation



The shoulder joint sling OmoSupport.easy is simple to apply and excellently adjustable to the individual body shape of the patient, due to its variable construction. The shoulder joint sling OmoSupport.easy is applied for injuries of the shoulder and immobilises the entire arm in a defined inner rotation position. By means of the hook- and loop closures, the position is precisely adjustable, while the open construction provides easy wound inspection. The innovative hand pocket is simple to open by hook- and loop closures and can be closed with one hand

Indication for use

Post traumatic, post surgery or conservative immobilisation of the shoulder joint for distortions, contusions, luxations and fractures

Special features

- Simple size options, „2 fit for all“, therefore low stock keeping
- Simple handling; easy application of the shoulder sling without help
- Partially flexible body strap; elastic insert of the body strap provides unrestricted breathing
- Innovative hand pocket closure; individual adjustment even for bandaged limbs
- Upper material with clinging structure for hook straps; individual adjustment

Item-No.	Size	Hip circumference	Colour
105R1/M	M	up to 110 cm	silver grey/black
105R1/L	L	110 – 140 cm	silver grey/black

Further shoulder joint orthoses:

OmoSupport.gilchrist (110R1/..), OmoSupport.air (100R1)

